



## SARS-CoV-2 Antigen Rapid Test

**Offers the Best Protection for  
Your Staff and Your Enterprise**

**Can accurately detect Omicron and Delta viruses**

- ✓ Reasonably Priced ✓ High Accuracy
- ✓ Approved by European Commission ✓ CE Certified
- ✓ Nasal Sample Collection ✓ Sample Disposal Bag Included

Germany's Most Authoritative  
Testing Institute  
Paul Ehrlich Institute  
Latest Test Result

$Cq \leq 25$  Sensitivity 100%  
 $Cq 25-30$  Sensitivity 95.7%

**Scored Top Spots in Overall Sensitivity  
Outstanding Results for Accuracy  
Widely Used by German Enterprises**

Hong Kong Distribution Agent:  
Prunus Jade Creation Limited (Established in 2016)  
[www.jademedicals.com](http://www.jademedicals.com)  
Telephone: 5508-1331



# 新冠病毒快速測試

給您的員工和企業最好的保障

能準確檢測Omicron及Delta病毒

- ✓ 價錢合理
- ✓ 準確度高
- ✓ 歐盟認可
- ✓ CE
- ✓ 鼻拭子取樣
- ✓ 內備樣品棄置袋方便安全

德國最權威檢測所  
Paul Ehrlich Institute  
最新測驗結果

Cq ≤ 25 敏感度 100%  
Cq 25-30 敏感度 95.7%

總敏感度前10名  
成績斐然  
德國巨企廣泛使用

香港代理:

Prunus Jade Creation Limited (成立於 2016 年)

[www.jademedicals.com](http://www.jademedicals.com)

電話: 5508-1331





# CERTIFICATE

**EC Certificate No. 1434-IVDD-526/2021**

EC Design-examination  
Directive 98/79/EC concerning  
*in vitro* diagnostic medical devices

Polish Centre for Testing and Certification certifies  
that manufactured by:

**Shenzhen Watmind Medical Co., Ltd**  
**Kengzi Subdistrict, Pingshan District No.16-1, 8th Floor,**  
**Building A, No.16-1, Jinhui Road, Jinsha Community,**  
**Kengzi Subdistrict, Pingshan District, 518118, Shenzhen, China**

in vitro diagnostic medical devices  
for self-testing

*The list of medical devices covered by this certificate is provided in the Annex 1*

In terms of design documentation, comply with requirements  
of Annex III (Section 6) to Directive 98/79/EC (as amended)  
Implemented into Polish law,  
as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 29.12.2021 to 27.05.2025

The date of issue of the Certificate: 29.12.2021

The date of the first issue of the Certificate: 29.12.2021

**CE** 1434

Issued under the Contract No. MD-67/2021  
Application No: 135/2021  
Certificate bears the qualified signature:  
Warsaw, 29/12/2021  
Module A1  
FBM-30-E\_10

President  
  
Aleksandra Kostrzewa

Digitally signed  
by Aleksandra  
Kostrzewa

President

**SARS-CoV-2 Ag Self-test Kit (Nasal Swab)**  
Instructions for Use (IFU)

www.mind

English

REF: LFA0401-LAE, LFA0401-GAR, LFA0401-ZAE

CE 1434

For lay use  
For in vitro diagnostic use only  
For use with nasal swab

This instruction for use (IFU) must be read carefully prior to use. Instruction for use must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions for use.

**INTENDED USE**

This kit is used for in vitro qualitative detection of Nucleocapsid (N) Protein antigen from SARS-CoV-2 in human nasal swab samples. This kit is authorized for lay use with self-collected observed direct nasal swab samples from individuals aged 15 years or older or adult collected nasal swab samples from individuals aged four years or older. Results are for the identification of SARS-CoV-2 Nucleocapsid Protein antigen. Antigen is generally detectable in nasal swab during the acute phase of infection. Positive results indicate the presence of viral antigen, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

Negative results should be treated as presumptive and confirmation with a molecular assay. Negative results cannot exclude SARS-CoV-2 infection and should not be used as the sole basis for patient management decisions. Negative results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19.

Individuals who test negative and continue to experience COVID-19-like symptoms should seek follow up care from their healthcare provider.

**SUMMARY AND EXPLANATION OF THE TEST**

Coronaviruses are a large family of viruses which may cause illness in animals or humans. SARS-CoV-2 is an enveloped, single-stranded RNA virus of the 2 gene. The virus can cause mild to severe respiratory illness and has spread globally.

The SARS-CoV-2 Ag Test Kit is a rapid lateral flow immunoassay for the qualitative detection of SARS-CoV-2 directly from nasal swab, without viral transport media. The kit contains all components required to carry out an assay for SARS-CoV-2.

**PRINCIPLE OF THE PROCEDURE**

The SARS-CoV-2 Ag Test Kit is a lateral flow immunoassay for the qualitative determination of Nucleocapsid Protein of SARS-CoV-2 virus in human nasal swab samples.

SARS-CoV-2 antibody is immobilized in the test region on nitrocellulose membrane. If the specimen contains SARS-CoV-2 antigen, during the assay specimen is allowed to react with the colored conjugate (SARS-CoV-2 antibody-colloidal gold conjugate); the molecules then migrate chromatographically on the membrane by the capillary action. An SARS-CoV-2 positive specimen produces a distinct color band in the test region, formed by the specific antibody-antigen colored conjugate complex [Au-SARS-CoV-2-Ab]-(SARS-CoV-2-Ag)-(SARS-CoV-2-Ab)\*. Absence of this colored band in the test region suggests a negative result. A colored band always appears in the control region serving as procedural control regardless of the specimen contains SARS-CoV-2 or not.

**REAGENTS AND MATERIALS**

Materials provided

Specificity	LFA0401-LAE (1 Function only) (unit)	LFA0401-GAR (2 Functions components)	LFA0401-ZAE (2 Functions components)
Test card with membrane in a sealed foil pouch	1	1	25
Sample extraction solution	1	1	25
Disposable Nasal Swab (N) swab	1	1	25
Diagram	1	1	25
Tube	1	1	25
Individual swab bag	1	1	25
Instruction for use	1	1	1

Materials required but not provided:  
\* Clock, timer or stopwatch

**PRECAUTIONS**

- For in vitro diagnostic use.
- This product has been authorized only for the detection of nucleocapsid protein from SARS-CoV-2, not for any other viruses or pathogens.
- Proper sample collection and handling are essential for correct results.
- Do not touch swab tip when handling the swab sample.
- Leave test card sealed in foil pouch until just before use. Do not use if pouch is damaged or open.
- Do not use kit past its expiration date.
- Do not mix test card and sample extraction solution from different kit lots.
- All kit components are single use items. Do not use with multiple specimens. Do not reuse the used test card.

**STORAGE AND STABILITY**

Kits should be stored in 2°C-8°C, valid for 18 months, forbidden to store under 2°C and avoid using expired products. Manufacture date (MFD) and Expiry date (EXP) marked on the label.

**TEST PROCEDURE**

Do not open the pouch until you are ready to perform a test, and the single-use test is suggested to be used under low environmental humidity (RH<70%) within 1 hour.

Before performing the test, you must read the instruction manual of this product completely, and please balance the test card and sample extraction solution to room temperature (15°C-25°C) before the test. Do not perform the test until the reagent was equilibrated to room temperature (15°C-25°C) so as to avoid affecting the accuracy of the experimental results.

**Wash Your Hands before Testing**

Before you start testing, wash your hands or use hand sanitizer. Make sure your hands are dry before starting.



**1. Open your test kit and you should have:**

REF	1	2	3	4	5
LFA0401-LAE	1	1	1	1	1
LFA0401-GAR	5	5	5	5	5
LFA0401-ZAE	25	25	25	25	25

**2. Open Pouch and place the card on a clean, dry, flat surface.**  
NOTE: Do not touch any parts on inside of the card.



**1. Open the sample extraction solution sealed vial.**



**4. Drop all the liquid into the tube.**



**5. Open Swab**  
Open swab package and take the swab out.  
NOTE: Keep fingers away from swab end.



**6. Sample Collection Process.**

Do not use nasal spray within 10 minutes before collecting nasal swab. Carefully insert the entire absorbent tip of the swab into the nostril (about 1.5-2 cm). Gently sample the nasal wall by rotating the swab in a circular path five times against the nasal wall. Slowly remove swab from the nostril. (This step should take approximately 10 seconds, ensuring to collect mucus and cells). Repeat the same process with the same swab in the other nostril!

Note: Simply rubbing the swab against one part of the inside of the nose or leaving the swab in the nose for less than 15 seconds is not a proper technique and may result in an insufficient sample.



**7. Elution of samples from swab.**

Place the swab into the sample tube and then completely immerse the swab head in the sample. Vigorously mix the solution by rotating the swab horizontally against the side of the tube at least 10 times (while submerged) and squeeze the tube 5 times by hand to ensure that the sample on the sampling swab is fully eluted into the sample extraction buffer.



**8. Squeeze the swab head along the inner wall of the tube to keep the liquid in the tube as much as possible.**



**9. Discard the swab, cover the drip head and wrangle tube for 5-6 times.**



**10. Dispense 100µL (4drops) of the specimen into the sample well ("S" well) on the card.**



**11. Wait 15 minutes.**  
NOTE: Do not disturb card during this time. False results can occur if the card is disturbed/moved or test results are read before 15 minutes.



**12. Interpret the test results at 15-20 minutes.**  
Do not interpret the results after 20 minutes.

**INTERPRETATION OF TEST RESULTS**

There are three types of results possible:

- Positive**  
Both red/purple test band (T) and red/purple control band (C) appear in window.



Note: The red/purple band in the test area (T) can show the color depth. However, within the specified observation time, regardless of the color of the ribbon, even a very weak ribbon should be judged as a positive result.

**In case of a positive test result:**

- COVID-19 infection is currently suspected.
- Immediately contact physician/family physician or local health department.
- Follow local guidelines for self-isolation.
- Have a PCR confirmatory test performed.

**2. Negative**

Only the red/purple control band (C) appears in window. The absence of a test band (T) indicates a negative result.



**In case of a negative test result:**

- Continue to follow all applicable rules regarding contact with others and protective measures.
- Even if the test is negative, an infection may be present.
- In case of suspicion, reperform test after 1 - 2 days, as the coronavirus cannot be accurately detected at all stages of an infection.

### 3. Invalid

There should always be a red/purple control band (C) in the control region regardless of test result. If control band (C) is not seen, it indicates that the incorrect operation process or the kit has deteriorated or damaged.



In case of an invalid test result:

- Possibly caused by incorrect test performance
- Repeat the test
- If test results remain invalid, contact a physician or COVID-19 testing center

### Handling instructions/Actions after the test result

- The following reasons may cause false negative results:
  - Inappropriate sample collection, using other non-matching solution, sample transfer time is too long (more than half an hour), the volume of solution added when diluted the swab are too much, non-standardized dilution operation, low virus titer in the sample, these may all lead to false negative results
  - Disturbance in final process may lead to changes in antigen epitopes, leading to false negative results.
- Analyze the possibility of false positive results:
  - Inappropriate sample collection, using other non-matching solutions, non-standardized dilution operation, these may all lead to false positive results
  - Cross-contamination of samples may lead to false positive results
  - False negative result from nucleic acid
  - Analyze the possibility of invalid result
    - If the sample volume is not enough, the chromatography cannot be carried out successfully
    - The test card would invalid if the package was broken. The packaging status must be carefully checked before use.

### DISPOSE IN TRASH

After test is completed, all used test components should put in biohazard waste bag and seal it. The trash must be disposed of in your household waste or in accordance with local disposal regulations.

### Wash Your Hands after Testing

After completing all test steps, wash your hands or use hand sanitizer.



### LIMITATIONS

- The result of the product should not be taken as a confirmed diagnosis, for clinical reference only. Judgments should be made along with RT-PCR results, clinical symptoms, epidemic condition and further clinical data.
- Due to the limitation of the detection method, the negative result cannot exclude the possibility of infection. The positive result should not be taken as a confirmed diagnosis. Judgments should be made along with clinical symptoms and further diagnosis methods.
- Positive test results do not rule out co-infections with other pathogens.
- False negative results are more likely after 5 days or more of symptoms.
- Negative results, from patients with symptoms onset beyond 7 days, should be treated as presymptomatic and confirmations with a molecular assay, if necessary, for patient management, may be performed.
- This reagent can only qualitatively detect SARS-CoV-2 antigens in human nasal swab samples. It cannot determine the specific antigen content in the sample.
- The accuracy of the test depends on the sample collection process. Inappropriate sample collection will affect the test results.

1. False negative results may occur if swabs are stored in their paper sheath after specimen collection.

2. Negative test results are not intended to rule in other non-SARS viral or bacterial infections.

10. A negative test result may occur if the level of antigen in a sample is below the detection limit of the test.

11. Cross reactions may be occur due to the N protein in SARS has a high homology with the non-coronavirus (SARS-CoV-2). However, the interpretation of the results is not affected during seasons without SARS infection.

### PERFORMANCE CHARACTERISTIC

#### 1. Analytical Performance

##### 1.1. Limit of detection

This kit was confirmed to detect 1.0-100 TCID50/mL of SARS-CoV-2 which was isolated from USA-WA1/2020, Gamma-irradiated.

##### 1.2. Cross reactivity

The following viruses and other Microorganisms have no effect on the test results.

Parental Gene-Sequence	Test Concentration	Test Result
<b>Respiratory Syncytial Virus A</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Respiratory Syncytial Virus B</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Influenza Virus</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Adenovirus Type 7</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Adenovirus Type 7</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Rhino virus</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Varicella-zoster virus</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Human coronavirus OC-43</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Human coronavirus HKU1</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Human coronavirus NL63</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Rubivirus</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Influenza H</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Influenza A</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Metapneumovirus</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Equine Herpes Virus</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>MERS-CoV</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>HCoV-229E</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Coronavirus NL63</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Epstein-Barr virus</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction
<b>Other Microorganisms</b>	1.0 x 10 <sup>7</sup> TCID <sub>50</sub> /mL	No cross reaction

An in-silico analysis was performed using the Basic Local Alignment Search Tool (BLAST) managed by the National Center for Biotechnology Information (NCBI) for SARS-coronavirus SARS-coronavirus shows 90-92% homology across 100% of the nucleocapsid sequence. Therefore, cross-reactivity is highly likely.

#### 1.3. Interfering Substances

The following interfering substances have no effect on the test results.

Substance	Active Ingredient	Concentration	Test Result
<b>Antibiotics</b>	Streptomycin	1.0 mg/mL	No interference
<b>Small Molecules</b>	Glucose	10 mg/mL	No interference
<b>Non-Threat Based Agents</b>	Phenol	10 mg/mL	No interference
<b>Other Substances</b>	Hydrogen Peroxide, Methyl	0.1% w/v	No interference
<b>Antiviral Drug</b>	Remdesivir	10 mg/mL	No interference
<b>Antifungal Agents</b>	Fluconazole	1 mg/mL	No interference

#### 1.4. High Dose Hook Effect

No high dose hook effect was observed when tested with up to a concentration of 1.0 x 10<sup>8</sup> TCID<sub>50</sub>/mL of heat inactivated SARS-CoV-2 virus.

### 2. Clinical study

The clinical evaluation was performed to compare the results obtained with the SARS-CoV-2 Ag Test Kit and a comparative reference manufacturer polymerase chain reaction test (Novel Coronavirus(2019-nCoV) Nucleic Acid Diagnostic Kit/PCR-Fluorescence Probing) manufactured by Sansiro Biotech Inc). Among patients, there are 103 positive and 302 negative samples by RT-PCR confirmed. The presentation of the results of the SARS-CoV-2 Ag Test Kit is as follows:

CT value	Number of sample	2019-nCoV RT-PCR Results	SARS-CoV-2 antigen result (compared with RT-PCR)
<18	66	positive	100% (66/66) (95% CI: 94.8% to 100%)
18-24	173	positive	100% (173/173) (95% CI: 99.7% to 100%)
>24	305	negative	100% (305/305) (95% CI: 99.7% to 100%)

Days	Number of sample	2019-nCoV RT-PCR Results	SARS-CoV-2 antigen result (compared with RT-PCR)
<7	100	positive	92.0% (92/100) (95% CI: 87.4% to 96.6%)
7-14	173	positive	100% (173/173) (95% CI: 99.7% to 100%)
>14	33	positive	100% (33/33) (95% CI: 94.8% to 100%)

Sensitivity: 92.0% 95%CI: 87.4%-96.6% for CT values <16.  
Sensitivity: 91.51% 95%CI: 86.6%-96.47% for onset of symptoms within 7 days.  
Specificity: 99.00% 95%CI: 98.51%-99.50%

### 3. Human Factors Study

Walmed conducted a human factor's study to evaluate whether home user patients or caregivers (lay user) could perform the test and accurately interpret test results from the SARS-CoV-2 Ag Card.

In this study, a total of 50 lay users, age 18 and older with either good or corrected vision (far-sighted or near-sighted) participated in a 30-minute session including an introduction, a product overview, and simulated use cases of SARS-CoV-2 Ag Test Kit result interpretation. Participants were asked to read and interpret a panel of 7 different SARS-CoV-2 Ag Card test results, including high positive, low positive, negative and invalid.

60-50 participants described the process of reading and interpreting the test card results as being easy. However, 6-50 of the participants commented that it was difficult to see some of the finer line conditions. A total of 340 trials were recorded in this study. Participants were able to perceive and interpret the results correctly for 327 trials, or 96.4% of the time. Positive results with stronger intensity lines were easier to read than the positive lines with less intensity. After the human factor evaluation, participants were asked for their overall impressions of the instructional materials they were provided. Nearly all participants (49-50) thought the instructions were straightforward and easy to understand and follow.

### 4. Usability Study

Walmed conducted a study to evaluate whether a home user can read the instructions and successfully perform the test steps for the SARS-CoV-2 Ag Card test, including swab collection at home, and correctly interpreting the results.

200 home users, including individuals (n=100) and caregivers (n=100), participated in the study. Each individual or caregiver pair participated in a 30-minute session with an instruction. The usability evaluation session included one simulated use of the SARS-CoV-2 Ag Test Kit. 95.0% (196 out of 200) home users produced a valid result (all negative) and 4 participants produced an invalid result. (The causes of the invalid result were insufficient amount of reagent added, and damage to the test strip). 199 out of 200 participants interpreted their test result correctly and 4 participants interpreted their result incorrectly (while they perceived a faint line in the sample window (as positive) when there was none (all results were verified by the study moderator).

The individual home use group completed 95.0% (196/200) of the total tasks steps correctly. The caregiver home user group completed 97.4% (124/127) of the total tasks steps correctly. The most common use errors observed during critical tasks included incorrectly swab collection and contacting the test strip with the hands or with the surface. 95.5% (197 out of 200) of the home (individual and caregiver) participants had positive impressions of the SARS-CoV-2 Ag Test Kit. The test was perceived as being easy to use.

### SYMBOLS

	For In Vitro Diagnostic Use		This symbol indicates the product's catalog number
	This symbol indicates that the product should be kept away from sunlight		This symbol indicates the product's Manufacturer
	This symbol indicates the product's Manufacturing Date		This symbol indicates that the product should be stored between 2-8°C
	This symbol indicates the product's Expiry Date		This symbol indicates that the product is for single use only. Do not re-use.
	This symbol indicates that you should consult the instructions for use.		This symbol indicates the product's Batch Number
	This symbol indicates that the product should be kept dry		EU Authorized Representative
	CE Mark		

### SARS-CoV-2 Ag Self-test Kit (Nasal Swab)

**Shenzhen Walmed Medical Co., Ltd.**  
 5th Floor, Building 4, No.10-1, Jianshe Road, Jianshe Community, Xixiangdi Street, Nanshan District, 518102 Shenzhen, China  
 Tel: +86-755-83070760  
 Email: gzh@walmed.com  
 Website: http://www.walmed.com

**Shenzhen International Medical Equip. Co., Ltd. (Europe)**  
 Address: Dillenschloß 30, 10237 Hamburg, Germany  
 Tel: +49-40-52101733  
 Email: shenzhen@imede.com

### Swab

**Jiangsu Hanting Medical Technology Co., Ltd.**  
 100#1, #17th Jiangsu Road, Tianzhu District, 210017 Jiangsu P.R. China

**Luzhou Laboratory Co., Ltd. (Europe)**  
 Dillenschloß 1, 10237, Hamburg, Germany

Date of Approval and Amendment of IFU: 2021-12-05



